



Run Info

Host Name	GXB01275 (localhost)
Position	X4
Experiment Name	kokia
Sample ID	Kk
Run ID	e2f8daf7-f0a0-4d10-83c3-68977f28ed98
Acquisition ID(s)	bbe2887a8ccfefe507bcb82f38d3bd8878ee7885, 210e1958f245a6d04e1ffe816e414c7bd05acd06
Flow Cell Id	FAU07188
Start Time	June 10, 14:15
Run Length	3d 0h 3m

Run Summary

Reads Generated	1.7 M
Passed Bases	6.1 Gb
Failed Bases	2.13 Gb
Estimated Bases	9.56 Gb
Percentage Basecalled	102%


Run Parameters

Flow Cell Type	FLO-MIN112
Kit	SQK-LSK112
Initial bias voltage	-200 mV
FAST5 output	Enabled
FASTQ output	Enabled
BAM output	Disabled
Bulk file output	Disabled
Active channel selection	Enabled
Basecalling	Enabled
Specified run length	72 hours
FAST5 reads per file	4000
FAST5 output options	vbz_compress,fastq,raw
FASTQ reads per file	4000
FASTQ output options	compress
Mux scan period	1 hour 30 minutes
Reserved pores	0 %
Basecall model	dna_r10.4_e8.1_hac.cfg
Read filtering	min_qscore=9
Read splitting	enable=on


Versions

MinKNOW	21.11.7
MinKNOW Core	4.5.4
Bream	6.3.5
Guppy	5.1.13

Cumulative Output Reads




Cumulative Output Bases




Read Length Histogram Estimated Bases - Outliers Discarded

Estimated N50: 8.82 kb




Read Length Histogram Basecalled Bases - Outliers Discarded

Estimated N50: 7.76 kb




Read Length Histogram Estimated Bases

Estimated N50: 8.82 kb




Read Length Histogram Basecalled Bases


Estimated N50: 7.76 kb




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History




Translocation Speed



QScore



Disk Write Performance



Run Debug Messages

- Mux scan for flow cell FAU07188 has found a total of 173 pores. 149 pores available for immediate sequencing June 13, 13:47
- Performing Mux Scan June 13, 13:44
- Mux scan for flow cell FAU07188 has found a total of 184 pores. 157 pores available for immediate sequencing June 13, 12:14
- Performing Mux Scan June 13, 12:12
- Mux scan for flow cell FAU07188 has found a total of 196 pores. 163 pores available for immediate sequencing June 13, 10:42
- Performing Mux Scan June 13, 10:39
- Mux scan for flow cell FAU07188 has found a total of 188 pores. 164 pores available for immediate sequencing June 13, 09:09
- Performing Mux Scan June 13, 09:07
- Mux scan for flow cell FAU07188 has found a total of 208 pores. 175 pores available for immediate sequencing June 13, 07:37
- Performing Mux Scan June 13, 07:34
- Mux scan for flow cell FAU07188 has found a total of 197 pores. 166 pores available for immediate sequencing June 13, 06:04
- Performing Mux Scan June 13, 06:02
- Mux scan for flow cell FAU07188 has found a total of 210 pores. 180 pores available for immediate sequencing June 13, 04:32
- Performing Mux Scan June 13, 04:29
- Mux scan for flow cell FAU07188 has found a total of 219 pores. 186 pores available for immediate sequencing June 13, 02:59
- Performing Mux Scan June 13, 02:57
- Mux scan for flow cell FAU07188 has found a total of 219 pores. 182 pores available for immediate sequencing June 13, 01:27
- Performing Mux Scan June 13, 01:24
- Mux scan for flow cell FAU07188 has found a total of 263 pores. 205 pores available for immediate sequencing June 12, 23:54
- Performing Mux Scan June 12, 23:52
- Mux scan for flow cell FAU07188 has found a total of 241 pores. 199 pores available for immediate sequencing June 12, 22:22
- Performing Mux Scan June 12, 22:19
- Mux scan for flow cell FAU07188 has found a total of 292 pores. 233 pores available for immediate sequencing June 12, 20:49
- Performing Mux Scan June 12, 20:47
- Mux scan for flow cell FAU07188 has found a total of 299 pores. 235 pores available for immediate sequencing June 12, 19:17
- Performing Mux Scan June 12, 19:14
- Mux scan for flow cell FAU07188 has found a total of 245 pores. 198 pores available for immediate sequencing June 12, 17:44
- Performing Mux Scan June 12, 17:42
- Mux scan for flow cell FAU07188 has found a total of 235 pores. 184 pores available for immediate sequencing June 12, 16:12
- Performing Mux Scan June 12, 16:09
- Mux scan for flow cell FAU07188 has found a total of 245 pores. 195 pores available for immediate sequencing June 12, 14:39
- Performing Mux Scan June 12, 14:36
- Mux scan for flow cell FAU07188 has found a total of 292 pores. 218 pores available for immediate sequencing June 12, 13:06

- Performing Mux Scan June 12, 13:03
- Mux scan for flow cell FAU07188 has found a total of 287 pores. 214 pores available for immediate sequencing June 12, 11:32
- Performing Mux Scan June 12, 11:30
- Mux scan for flow cell FAU07188 has found a total of 302 pores. 217 pores available for immediate sequencing June 12, 09:59
- Performing Mux Scan June 12, 09:56
- Mux scan for flow cell FAU07188 has found a total of 304 pores. 217 pores available for immediate sequencing June 12, 08:25
- Performing Mux Scan June 12, 08:23
- Mux scan for flow cell FAU07188 has found a total of 324 pores. 221 pores available for immediate sequencing June 12, 06:52
- Performing Mux Scan June 12, 06:49
- Mux scan for flow cell FAU07188 has found a total of 327 pores. 225 pores available for immediate sequencing June 12, 05:18
- Performing Mux Scan June 12, 05:16
- Mux scan for flow cell FAU07188 has found a total of 341 pores. 220 pores available for immediate sequencing June 12, 03:45
- Performing Mux Scan June 12, 03:42
- Mux scan for flow cell FAU07188 has found a total of 354 pores. 221 pores available for immediate sequencing June 12, 02:11
- Performing Mux Scan June 12, 02:09
- Mux scan for flow cell FAU07188 has found a total of 372 pores. 240 pores available for immediate sequencing June 12, 00:38
- Performing Mux Scan June 12, 00:35
- Mux scan for flow cell FAU07188 has found a total of 387 pores. 252 pores available for immediate sequencing June 11, 23:04
- Performing Mux Scan June 11, 23:02
- Mux scan for flow cell FAU07188 has found a total of 428 pores. 277 pores available for immediate sequencing June 11, 21:31
- Performing Mux Scan June 11, 21:28
- Mux scan for flow cell FAU07188 has found a total of 420 pores. 264 pores available for immediate sequencing June 11, 19:57
- Performing Mux Scan June 11, 19:55
- Mux scan for flow cell FAU07188 has found a total of 450 pores. 274 pores available for immediate sequencing June 11, 18:24
- Performing Mux Scan June 11, 18:21
- Mux scan for flow cell FAU07188 has found a total of 499 pores. 303 pores available for immediate sequencing June 11, 16:50
- Performing Mux Scan June 11, 16:48
- Mux scan for flow cell FAU07188 has found a total of 537 pores. 319 pores available for immediate sequencing June 11, 15:17
- Performing Mux Scan June 11, 15:14
- Mux scan for flow cell FAU07188 has found a total of 562 pores. 342 pores available for immediate sequencing June 11, 13:43
- Performing Mux Scan June 11, 13:41
- Mux scan for flow cell FAU07188 has found a total of 554 pores. 339 pores available for immediate sequencing June 11, 12:10
- Performing Mux Scan June 11, 12:08
- Mux scan for flow cell FAU07188 has found a total of 620 pores. 356 pores available for immediate sequencing June 11, 10:36
- Performing Mux Scan June 11, 10:34

- Mux scan for flow cell FAU07188 has found a total of 651 pores. 372 pores available for immediate sequencing June 11, 09:03
- Performing Mux Scan June 11, 09:01
- Mux scan for flow cell FAU07188 has found a total of 688 pores. 390 pores available for immediate sequencing June 11, 07:30
- Performing Mux Scan June 11, 07:27
- Mux scan for flow cell FAU07188 has found a total of 726 pores. 406 pores available for immediate sequencing June 11, 05:56
- Performing Mux Scan June 11, 05:54
- Mux scan for flow cell FAU07188 has found a total of 752 pores. 406 pores available for immediate sequencing June 11, 04:23
- Performing Mux Scan June 11, 04:20
- Mux scan for flow cell FAU07188 has found a total of 767 pores. 405 pores available for immediate sequencing June 11, 02:49
- Performing Mux Scan June 11, 02:47
- Mux scan for flow cell FAU07188 has found a total of 836 pores. 413 pores available for immediate sequencing June 11, 01:16
- Performing Mux Scan June 11, 01:13
- Mux scan for flow cell FAU07188 has found a total of 848 pores. 427 pores available for immediate sequencing June 10, 23:42
- Performing Mux Scan June 10, 23:40
- Mux scan for flow cell FAU07188 has found a total of 876 pores. 435 pores available for immediate sequencing June 10, 22:09
- Performing Mux Scan June 10, 22:06
- Mux scan for flow cell FAU07188 has found a total of 927 pores. 442 pores available for immediate sequencing June 10, 20:35
- Performing Mux Scan June 10, 20:33
- Mux scan for flow cell FAU07188 has found a total of 971 pores. 461 pores available for immediate sequencing June 10, 19:02
- Performing Mux Scan June 10, 18:59
- Mux scan for flow cell FAU07188 has found a total of 996 pores. 465 pores available for immediate sequencing June 10, 17:28
- Performing Mux Scan June 10, 17:26
- Mux scan for flow cell FAU07188 has found a total of 1060 pores. 483 pores available for immediate sequencing June 10, 15:55
- Performing Mux Scan June 10, 15:52
- Mux scan for flow cell FAU07188 has found a total of 1084 pores. 491 pores available for immediate sequencing June 10, 14:21
- Performing Mux Scan June 10, 14:19
- Starting sequencing procedure June 10, 14:19
- Waiting up to 300 seconds for temperature to stabilise at 30.0°C June 10, 14:15