



## Run Info

Host Name	GXB01275 (localhost)
Position	X1
Experiment Name	kokia
Sample ID	Kc
Run ID	325fe942-7810-47ec-815b-accb048c1ff4
Acquisition ID(s)	6dcbb864a508f1e9d4b570e5265cebbba65ca4dad, d8294921947fa9e9a8afb74474e7dda80aafc71d
Flow Cell Id	FAU08661
Start Time	June 10, 14:15
Run Length	3d 0h 3m

## Run Summary

Reads Generated	1.89 M
Passed Bases	9.01 Gb
Failed Bases	2.65 Gb
Estimated Bases	13.4 Gb
Percentage Basecalled	102%


## Run Parameters

Flow Cell Type	FLO-MIN112
Kit	SQK-LSK112
Initial bias voltage	-200 mV
FAST5 output	Enabled
FASTQ output	Enabled
BAM output	Disabled
Bulk file output	Disabled
Active channel selection	Enabled
Basecalling	Enabled
Specified run length	72 hours
FAST5 reads per file	4000
FAST5 output options	vbz_compress,fastq,raw
FASTQ reads per file	4000
FASTQ output options	compress
Mux scan period	1 hour 30 minutes
Reserved pores	0 %
Basecall model	dna_r10.4_e8.1_hac.cfg
Read filtering	min_qscore=9
Read splitting	enable=on


## Versions

MinKNOW	21.11.7
MinKNOW Core	4.5.4
Bream	6.3.5
Guppy	5.1.13

Cumulative Output Reads




Cumulative Output Bases




**Read Length Histogram Estimated Bases - Outliers Discarded**

Estimated N50: 10.71 kb




**Read Length Histogram Basecalled Bases - Outliers Discarded**

Estimated N50: 9.34 kb




**Read Length Histogram Estimated Bases**

Estimated N50: 10.71 kb




**Read Length Histogram Basecalled Bases**


Estimated N50: 9.34 kb




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History



Translocation Speed




QScore





Disk Write Performance



**Run Debug Messages**

- Mux scan for flow cell FAU08661 has found a total of 163 pores. 142 pores available for immediate sequencing June 13, 13:53
- Performing Mux Scan June 13, 13:51
- Mux scan for flow cell FAU08661 has found a total of 173 pores. 146 pores available for immediate sequencing June 13, 12:21
- Performing Mux Scan June 13, 12:19
- Mux scan for flow cell FAU08661 has found a total of 206 pores. 168 pores available for immediate sequencing June 13, 10:48
- Performing Mux Scan June 13, 10:46
- Mux scan for flow cell FAU08661 has found a total of 187 pores. 155 pores available for immediate sequencing June 13, 09:16
- Performing Mux Scan June 13, 09:14
- Mux scan for flow cell FAU08661 has found a total of 224 pores. 187 pores available for immediate sequencing June 13, 07:43
- Performing Mux Scan June 13, 07:41
- Mux scan for flow cell FAU08661 has found a total of 221 pores. 185 pores available for immediate sequencing June 13, 06:11
- Performing Mux Scan June 13, 06:09
- Mux scan for flow cell FAU08661 has found a total of 237 pores. 190 pores available for immediate sequencing June 13, 04:38
- Performing Mux Scan June 13, 04:36
- Mux scan for flow cell FAU08661 has found a total of 260 pores. 213 pores available for immediate sequencing June 13, 03:06
- Performing Mux Scan June 13, 03:04
- Mux scan for flow cell FAU08661 has found a total of 254 pores. 192 pores available for immediate sequencing June 13, 01:33
- Performing Mux Scan June 13, 01:31
- Mux scan for flow cell FAU08661 has found a total of 283 pores. 221 pores available for immediate sequencing June 13, 00:00
- Performing Mux Scan June 12, 23:57
- Mux scan for flow cell FAU08661 has found a total of 313 pores. 234 pores available for immediate sequencing June 12, 22:27
- Performing Mux Scan June 12, 22:24
- Mux scan for flow cell FAU08661 has found a total of 338 pores. 247 pores available for immediate sequencing June 12, 20:53
- Performing Mux Scan June 12, 20:51
- Mux scan for flow cell FAU08661 has found a total of 311 pores. 223 pores available for immediate sequencing June 12, 19:20
- Performing Mux Scan June 12, 19:17
- Mux scan for flow cell FAU08661 has found a total of 350 pores. 239 pores available for immediate sequencing June 12, 17:46
- Performing Mux Scan June 12, 17:44
- Mux scan for flow cell FAU08661 has found a total of 303 pores. 211 pores available for immediate sequencing June 12, 16:13
- Performing Mux Scan June 12, 16:10
- Mux scan for flow cell FAU08661 has found a total of 332 pores. 223 pores available for immediate sequencing June 12, 14:39
- Performing Mux Scan June 12, 14:37
- Mux scan for flow cell FAU08661 has found a total of 363 pores. 249 pores available for immediate sequencing June 12, 13:06

- Performing Mux Scan June 12, 13:03
- Mux scan for flow cell FAU08661 has found a total of 397 pores. 257 pores available for immediate sequencing June 12, 11:32
- Performing Mux Scan June 12, 11:30
- Mux scan for flow cell FAU08661 has found a total of 444 pores. 281 pores available for immediate sequencing June 12, 09:59
- Performing Mux Scan June 12, 09:56
- Mux scan for flow cell FAU08661 has found a total of 477 pores. 312 pores available for immediate sequencing June 12, 08:25
- Performing Mux Scan June 12, 08:23
- Mux scan for flow cell FAU08661 has found a total of 499 pores. 313 pores available for immediate sequencing June 12, 06:52
- Performing Mux Scan June 12, 06:49
- Mux scan for flow cell FAU08661 has found a total of 461 pores. 285 pores available for immediate sequencing June 12, 05:18
- Performing Mux Scan June 12, 05:16
- Mux scan for flow cell FAU08661 has found a total of 496 pores. 299 pores available for immediate sequencing June 12, 03:45
- Performing Mux Scan June 12, 03:43
- Mux scan for flow cell FAU08661 has found a total of 552 pores. 329 pores available for immediate sequencing June 12, 02:11
- Performing Mux Scan June 12, 02:09
- Mux scan for flow cell FAU08661 has found a total of 585 pores. 350 pores available for immediate sequencing June 12, 00:38
- Performing Mux Scan June 12, 00:36
- Mux scan for flow cell FAU08661 has found a total of 580 pores. 351 pores available for immediate sequencing June 11, 23:04
- Performing Mux Scan June 11, 23:02
- Mux scan for flow cell FAU08661 has found a total of 620 pores. 351 pores available for immediate sequencing June 11, 21:31
- Performing Mux Scan June 11, 21:29
- Mux scan for flow cell FAU08661 has found a total of 669 pores. 369 pores available for immediate sequencing June 11, 19:57
- Performing Mux Scan June 11, 19:55
- Mux scan for flow cell FAU08661 has found a total of 655 pores. 361 pores available for immediate sequencing June 11, 18:24
- Performing Mux Scan June 11, 18:22
- Mux scan for flow cell FAU08661 has found a total of 682 pores. 370 pores available for immediate sequencing June 11, 16:50
- Performing Mux Scan June 11, 16:48
- Mux scan for flow cell FAU08661 has found a total of 725 pores. 400 pores available for immediate sequencing June 11, 15:17
- Performing Mux Scan June 11, 15:15
- Mux scan for flow cell FAU08661 has found a total of 775 pores. 402 pores available for immediate sequencing June 11, 13:44
- Performing Mux Scan June 11, 13:41
- Mux scan for flow cell FAU08661 has found a total of 795 pores. 400 pores available for immediate sequencing June 11, 12:10
- Performing Mux Scan June 11, 12:08
- Mux scan for flow cell FAU08661 has found a total of 806 pores. 407 pores available for immediate sequencing June 11, 10:37
- Performing Mux Scan June 11, 10:34

- Mux scan for flow cell FAU08661 has found a total of 858 pores. 423 pores available for immediate sequencing June 11, 09:03
- Performing Mux Scan June 11, 09:01
- Mux scan for flow cell FAU08661 has found a total of 892 pores. 433 pores available for immediate sequencing June 11, 07:30
- Performing Mux Scan June 11, 07:27
- Mux scan for flow cell FAU08661 has found a total of 946 pores. 442 pores available for immediate sequencing June 11, 05:56
- Performing Mux Scan June 11, 05:54
- Mux scan for flow cell FAU08661 has found a total of 982 pores. 459 pores available for immediate sequencing June 11, 04:23
- Performing Mux Scan June 11, 04:20
- Mux scan for flow cell FAU08661 has found a total of 1041 pores. 464 pores available for immediate sequencing June 11, 02:49
- Performing Mux Scan June 11, 02:47
- Mux scan for flow cell FAU08661 has found a total of 1095 pores. 473 pores available for immediate sequencing June 11, 01:16
- Performing Mux Scan June 11, 01:13
- Mux scan for flow cell FAU08661 has found a total of 1119 pores. 475 pores available for immediate sequencing June 10, 23:42
- Performing Mux Scan June 10, 23:40
- Mux scan for flow cell FAU08661 has found a total of 1172 pores. 482 pores available for immediate sequencing June 10, 22:09
- Performing Mux Scan June 10, 22:06
- Mux scan for flow cell FAU08661 has found a total of 1216 pores. 484 pores available for immediate sequencing June 10, 20:35
- Performing Mux Scan June 10, 20:33
- Mux scan for flow cell FAU08661 has found a total of 1257 pores. 489 pores available for immediate sequencing June 10, 19:02
- Performing Mux Scan June 10, 18:59
- Mux scan for flow cell FAU08661 has found a total of 1296 pores. 502 pores available for immediate sequencing June 10, 17:28
- Performing Mux Scan June 10, 17:26
- Mux scan for flow cell FAU08661 has found a total of 1352 pores. 499 pores available for immediate sequencing June 10, 15:55
- Performing Mux Scan June 10, 15:52
- Mux scan for flow cell FAU08661 has found a total of 1413 pores. 506 pores available for immediate sequencing June 10, 14:21
- Performing Mux Scan June 10, 14:19
- Starting sequencing procedure June 10, 14:19
- Waiting up to 300 seconds for temperature to stabilise at 30.0°C June 10, 14:15